




Factory door screening of food products: the potential and limitations of profiling techniques such as NIR

Dr Gerard Downey
Teagasc
Ashtown Food Research Centre, Dublin 15

Traceability System Overview



Factory Door Movements



What do we want to know about these materials *vis-à-vis* authenticity?

- **Conformance to specification** on purchase order or delivery documentation in the case of incoming material
- **Conformance to specification** on food label in the case of outgoing material

How do we address these knowledge requirements?

- Chemical or physical analysis
 - Time-consuming, expensive, destructive etc.



- Profiling methods
 - Generally fast, may be inexpensive, may be non-destructive, limited or no reagent disposal issues



Profiling Methods

- Vibrational spectroscopy
 - Near infrared, mid-infrared, Raman
- Separation methods
 - Combinations of LC, GC with mass spectroscopy (MS); IRMS
- Nuclear magnetic resonance (NMR)

TRACE: Tracing Food Commodities in Europe








www.trace.eu.org

trace.enquiries@trace.eu.org

Information Type and Specificity of Profiling Methods

Method	Information	Specificity
NIR	Molecular bond vibration – overtones + combinations	Low to medium
FT-IR	Molecular bond vibration – fundamental	Medium
Raman	Light scatter	Low to medium
*-MS	Molecular fragment mass	High
NMR	Chemical shift	High




Speed of Analysis and Ease of Industrial Deployment of Profiling Methods

Method	Speed	Ease of Industrial Deployment
NIR	Fast - seconds	Very easy
FT-IR	Fast – seconds (ATR)	Easy
Raman	Fast - seconds	Moderate
*-MS	Slow	Low
NMR	Slow	Low

What type of questions do we ask of the analytical data?

- Is this sample the same as some specified material?
 - **Classification**
- To which of a limited number of possible sample types does this sample belong?
 - **Discrimination**



Classification and Discrimination

- Need a **comprehensive** collection of **representative** samples of material to develop **robust and well-defined** models
- For naturally-occurring samples, **several harvest seasons** need to be included
- Whether samples for test are natural or processed, **sample collection** period can be **prolonged**

Classification and Discrimination

- Correct use requires significant knowledge of and familiarity with **chemometric procedures**
- How to **describe confidence** regarding outputs from predictions of real, unknown samples?

Classification Options

➤ SIMCA (soft independent modelling of class analogy)

- models a given type of *e.g.* raw material or finished product
- projects unknown samples into model space
- designates unknown as belonging or not belonging to model
- probability of designation is statistically derived

SIMCA Deployment

➤ Example of **honey adulteration**

- Irish artisanal honey (n=580)
- adulterated by fully-inverted beet syrup (n=280), high fructose corn syrup (n=160), partial invert cane syrup (n=120), dextrose syrup (n=160) and beet sucrose (n=120)
- all adjusted to 70 °Brix with distilled water; FT-IR ATR spectra collected
- SIMCA model (using 6 PCs) was developed on 290 authentic honeys; test sample file contained 290 authentic and 800 adulterated samples

SIMCA Honey Results

Sample type	No. of samples	No. correctly identified	No. incorrectly identified	% correctly identified	% incorrectly identified	
Honey	290	279	11	96.2	3.8	✓
BS Adulterated	120	117	3	97.5	2.5	✓
CS adulterated	120	115	5	95.8	4.2	✓

SIMCA Honey Results

Sample type	No. of samples	No. correctly Identified	No incorrectly Identified	% correctly Identified	% incorrectly identified	
Honey	290	279	11	96.2	3.8	
PICS Adulterated	120	62	58	51.2	48.3	✘
HFCS Adulterated	160	2	158	1.3	98.8	✘
IB adulterated	280	19	241	6.8	93.2	✘

Reasons for Failure

- Solutions are chemically identical, *or*
- Information on chemical differences is not contained in FT-IR ATR spectra, *or*
- Chemometric operations are insufficient to extract information

Discrimination

- Often stated as a limited class problem
- An example could be does this sample match one of a limited number of types?
- Some recent work involving NIR analysis of honey may serve as an example

Geographic origin of honey

- Preliminary feasibility study
- Sample set: 125 filtered honeys
 - 25 Irish
 - 25 Argentinean
 - 50 Czech
 - 25 Hungarian



Woodcock, T., Downey, G., Kelly, J.D.K. and O'Donnell, C.
J Ag. Fd Chem., 555(22), 9128-9134.

Four class PLS2 model results

Data	Argentinean (n=8)			Czech (n=16)		Hungarian (n=8)		Irish (n=8)	
	#L	% correct	% false	% correct	% false	% correct	% false	% correct	% false
Pretreatment		classification	positives	classification	positives	classification	positives	classification	positives
Raw Data	4	50	0	100	0	100	0	75	37.5
1st der, 5 pts	5	37.5	0	100	0	100	0	100	25
1st der, 9 pts	5	37.5	0	100	0	100	0	100	25
1st der, 13 pts	5	37.5	0	100	0	100	0	100	25
1st der, 21 pts	5	25	0	100	0	100	0	100	25
2nd der, 5 pts	5	62.5	0	100	0	100	0	87.5	50
2nd der, 9 pts	5	62.5	0	100	0	100	0	87.5	50
2nd der, 13 pts	5	50	0	100	0	100	0	100	37.5
2nd der, 21 pts	5	50	0	100	0	100	0	100	37.5
SNV	5	75	0	100	0	100	0	87.5	75


Argentine and Irish Separation

Data	#L	Argentinean (n=8)		Irish (n=8)	
		% correct	% false	% correct	% false
Pretreatment		classification	positives	classification	positives
Raw Data	5	87.5	0	100	12.5
1st der, 5 pts	4	100	0	100	0
1st der, 9 pts	4	100	0	100	0
1st der, 13 pts	4	100	0	100	0
1st der, 21 pts	4	100	0	100	0
2nd der, 5 pts	4	100	0	100	0
2nd der, 9 pts	4	100	0	100	0
2nd der, 13 pts	4	100	0	100	0
2nd der, 21 pts	4	100	0	100	0
SNV	4	100	0	100	0

Operation of PLS1


- Ascribe dummy variables to each sample type e.g.
 - 1 for honey belonging to a given class, &
 - 0 for other honeys
- Predict unknown honeys
- Select cut-off value for class ascription e.g. 0.5

Practical Outputs




Utilisation

- Such a model may be used to generate a **specification** derived from multivariate data
- How to decide on **value of cut-off** for specification?



How does the model work?


- Examination of e.g. PLS loading coefficients
- Example shown is from Argentine:Irish PLS1 model
- **Difficulty in unambiguous identification** of chemical species is apparent



Variable Selection

- Can often improve especially SIMCA models
- Permit selection of limited number of variables grouped in such a way as to provide a single number cut-off value
- Particularly applicable to $^1\text{H}/\text{MS}$ and NMR data

Variable Selection + Genetic Programming



Exploitation Scenario - 1


- Develop classification or discriminant model for material of interest using profiling technique
- Ensure robustness of model by creating extensive and representative spectral database
- Continually update as poorly predicted authentic samples are encountered

Exploitation Scenario - 2

- Each “producer” can generate fingerprint model of his/her product
- Production batches of this product may be checked using the model and guaranteed as meeting specification
- **Caveat:** In the absence of identifiable marker compounds such specifications may be used only to confirm conformance – they are not exclusive
- There is a trade-off between speed and accuracy; profiling methods can optimise this trade-off

Conclusions

- Fingerprinting and profiling methods can provide sufficient reassurance for efficient trading
- Industry needs to grasp the opportunity and invest in database collection by some technique
- Positive results from conformance tests may be used for marketing purposes as well as QC checks



C'est fini!

*Merci pour votre
attention.*